Mikroskope *Microscope*

2008/09

BRESSER®

Liebe Mikroskop-Interessentin, lieber Mikroskop-Interessent,

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Introduction

Dear customer,

you now have the new Bresser microscope catalogue in your hands. It contains many completely new microscopes. In order to make the right choice in choosing the correct microscope, please take a few minutes to read through our introductory text with explanations about the different types of microscope.

Advanced amateurs and professionals can refer directly to the technical data of the new Bresser Science line; these products are introduced commencing from page 24.

Generally, there are two main types of microscopes: Incident light and transmission type.

An **incident light microscope** is suitable for opaque samples like coins, stamps, minerals or insects. Here, the sample is illuminated from above and the reflected light is magnified by the objective. Stereo microscopes like the Biolux ICD or Researcher ICD give a three-dimensional view due to their twin objectives.

With **transmission-type microscopes**, the sample is penetrated by the light. In order to make that possible, the sample must be transparent, be prepared and put on a blank slide. Because such objects are mostly from within the realm of biology, these microscopes are also called "biological microscopes". Such thin sample slices give a fascinating view inside organs, cell cultures or even minerals.

Additionally, there are microscopes that are equipped with **both kinds of illumination**. Two examples are the Duolux and the Researcher ICD.

The Bresser microscope portfolio is divided in **three lines**: Models for beginners (starting at page 10), for school and hobbyists (from page 15 on) and advanced systems for science and research (from page 24 on).

The **beginner models** are very suitable for children and occasional observers. They provide an interesting entry into the "small" world.

The "school & hobby" line offers solid optical and mechanical quality at a reasonable price. The main emphasis is lies on practical use "always and everywhere". Many of these microscopes are equipped with internal rechargeable batteries for outdoor-use. This product line is ideal for schools and hobbyists. It is versatile and expandable with many useful accessories. These microscopes can therefore grow with the demands of the user.

The new "Science" line builds the base for demanding amateurs and professional observers like universities, scientific research and environment institutions. Precise high-class optics and mechanical components meet the high demands of daily laboratory use. Special illumination systems make specialist applications like fluorescence or polarisation microscopy possible or even inverse observation of sediments etc. A special accessory programme rounds out this product segment (see Microscope Matrix on pages 8-9).

Table of Contents:

Important points on microscopy	6-7
Microscope Matrix	8-9
Beginner Models	
Biotar	10
Microscope Set	11
Duolux	12-13
Biolux ICD / Halogen	14
Models for School and Hobby	
Biorit	15
Erudit MO	16
Erudit DLX	17
Researcher Bino + Trino	18-19
New ICD	20
Free Page	21
Researcher ICD	22-23
"Science" Models	
Science TRM-301	24-25
Science ADL-601P	26-27
Science ADL-601F	28-29
Science MPO-401	30-31
Science IVM-401	32-33
Science MTL-201	34-35
Science ETD-101	36-37
Accessories	
MicrOcular	38
Optional accessories	39-40
Glossary	
Microscope ABC's	47-49
Company profile	51

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Important points of microscopy

What characterises a modern microscope?

Two lenses - eyepiece and lens itself. If the lens is of poor quality then the eyepiece will magnify this lack as well. That's why good lenses are usually expensive.

Microscope Components

A modern microscope may consist of a variety of parts. See below. We have listed the main components in the following image:


Microscopy – what is it?

Micro comes from the Greek and means, „small“. Scopy is also of Greek origin and means, „look“.

How do I observe properly?

It's important the microscope be firmly fixed on a flat stable surface. The slide must be under the lens on the slide plate.


Always use the smallest eyepiece for initial adjustment, together with the weakest lens. Progress to coarse and fine focussing.

The accessories that are included with the microscope vary depending on their use, they may therefore be different than the ones depicted.

Lenses are built into microscopes – but how do they work?

The illustration shows how a lens works. The lens glass changes the direction of the light hitting the hammer. This makes it seem larger to the human eye (dotted lines).

Important points of microscopy

Further light bundling

The condenser (illustration below)

Bright lights generate a lot of heat that may damage specimens. To have bright light despite this a condenser can be used. The condenser consists of a lens system that bundles light onto the specimen to increase its brightness.

The diaphragm

- Iris diaphragm (illustration above)
Aperture size can be adjusted using the lateral lever as desired. It can be continuously adjusted from large to small.
- Aperture diaphragm (illustration below) A disc with holes of differing sizes can be used to adjust light quantity by turning it.

How to focus correctly

Adjust the distance between lens and specimen; always adjust from shortest to longest.

This means the lens is always to be made more distant from the slide.

Use coarse adjustment until an image appears. Then fine adjustment until the image is properly focussed.

Fine adjustment moves the lens in tiny parts of a mm to or from the specimen.

Note before buying

There are many different kinds of microscope. Before buying make sure you know the purpose you want your microscope to fulfil in the long term.

- Student microscope
Economical model with relatively low optical quality.
- Hobby microscope
Higher quality, good cost/benefit ratio. May accept some accessories to improve capacity.
- Professional microscope
Very high optical capacity. Can be improved by adding accessories.

Microscopes aren't judged by their magnification capacity alone. The quality of the lens and eyepiece/s is much more important.

Look for the lens apertures. They tell you what the lenses are individually capable of.

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

The Bresser beginner microscope BIOTAR has an excellent price-performance ratio. It includes numerous accessories, that make entry into the world of microscopy easy. The hard carry case offers perfect storage space for the accessories. To discover this „new“ world. The BIOTAR has a die-cast body and offers great optics.

Included with your purchase:

Hard carry case, Zoom eyepiece 10x - 20x, Objectives (5x, 20x, 60x), microtome, projector, hatchery tank, shrimp eggs, yeast (shrimp feed), sea salt, pipette, 10 blank slides, 5 prepared slides, 10 cover glasses, replacement bulb, projector, microscope utensils and much more ...

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High grade metal bodied student microscope with barlow lens zoom system, which magnifies the image up to an additional 1.6x and the Microcular I (PC eyepiece with 320x240 Pixels resolution).

Included with your purchase:

Barlow lens, Microcular I, Ulead Photo Explorer image editing software, USB connection cable, 2 optical eyepieces, AC power adapter, 5 preparation slides, 10 blank slides, 20 covering glasses and dust cap.

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Optional accessories can be found on pages 39 to 40.

Art.-Nr. · Art. No.:	50-11100
Modell · Model:	Mikroskop Set 40-1024x
Typ · Type:	Durchlicht · <i>Transmission</i>
Vergrößerung · <i>Magnification</i> :	40x - 1024x
Okulare · <i>Eyepieces</i> :	WF 10x / WF 16x
Objektive · <i>Objectives</i> :	4x / 10x / 40x
Beleuchtung · <i>Illumination</i> :	LED Licht · <i>LED light</i>

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The DuoLux is a complete new Bresser development. It combines the advantages of an incident light type with a transmission-type microscope to a handy combined microscope. With an integrated slide table, the classical microscope stage can be quickly transformed into an incident light stage with a 38mm diameter work field.

The 45mm variable distance between objective and work field is large enough to be able to view most interesting objects with incident light. The illumination source can easily be adjusted with a little joystick – also a new development by Bresser. Both illuminators can be switched on and off and are separately dimmable and can be used with battery or AC adaptor. For transmission-type microscopy, the integrated filter wheel is an optimal addition.

The complete standard accessories (e.g. electronic eyepiece, slides, preparation aids) can be stored in the integrated drawer.

Included with your purchase:

DIN Eyepieces (WF10x, WF16x), Objectives (4x, 10x, 40x), Microtome, Hatchery, Shrimp eggs, Yield (shrimp food), Sea salt, Pipette, 1.6x Barlow lens, hatchery, MicrOcular I, Ulead Photo Explorer image editing software (Windows 2000 oder XP), USB connection cable, AC power adapter*

- ❶ Wechselschieber · Slide table
- ❷ Zubehör-Schublade · Drawer for accessories
- ❸ MikrOkular I inklusive · MicrOcular I included
- ❹ Auswechselbare Filter/Linsen
Changeable filters/lenses
- ❺ Farbfilterscheibe · Colour filter pane

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This handy stereo microscope BIOLUX ICD with binocular head is ideal for mobile applications, due to its lightness. The illuminator is powered by batteries, which makes the microscope independent of any external electrical sources. The eye relief is adjustable, as with all bin- or trinocular Bresser microscopes. With optional available additional eyepieces, other magnifications than the standard 20x can be reached. The BIOLUX ICD, has a great optical quality, giving impressive images and optimal results in observations of coins, minerals, stamps, small electronic components, etc.*

Included with your purchase:

Eyepieces (WF10x), stereo objective 2x, box with 12 mineral samples.

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The Bresser BIORIT biological microscope is designed for the beginner in microscopy, for hobby and study. On the other hand, it has a superior specification. A part of this is the accu, which is mounted in the microscope base. The accu makes mobile use possible, as well as the AC adaptor which allows use in a single location.

You can observe the samples on your PC monitor with the Bresser MicrOcular I (320 x 240 pixels), which is included. Its software enables you to save and replay your microscope observations.*

Included with your purchase:

DIN eyepieces (WF10x, WF16x), Objectives (4x, 10x, 40x), mechanical desk with Nonius-scale, Barlow lens, MicrOcular I, Ulead Photo Explorer image editing software, USB connection cable, AC power adapter,

Art.-Nr. · Art. No.:	51-01000
Modell · Model:	Biorit
Typ · Type:	Durchlicht · Transmission
Vergrößerung · Magnification:	40x - 1280x
Okulare · Eyepieces:	DIN WF10x & WF16x (Paar · Pair)
Objektive · Objectives:	4x / 10x / 40x
Beleuchtung · Illumination:	LED Akku

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
The ERUDIT MO is a compact and versatile microscope for students and beginners. Equipped with accu-powered LED illumination, it is ready to go everywhere. The three objectives, in combination with three eyepieces and a Barlow-zoom-system, give a magnification range from 20x to 1536x.

The ERUDIT MO has a rigid metal body and comes with the accessories as standard equipment.

Included with your purchase:

Eyepieces (WF5x, WF10x, WF16x), objectives (4x, 10x, 16x), 1.6x Barlow lens, mechanical desk with Nonius-scale, Box with preparations, Microcular II, Ulead Photo Explorer image editing software, USB connection cable, wall adapter, alu travel case*

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The transmitted light illumination of the ERUDIT DLX is supplied by a LED with dimmer (for continuous variation of light intensity) and a height-adjustable condenser (with iris diaphragm and filter holder). This allows the instrument to have not only the 4x, 10x and 40x DIN objective, but also the 100x oil immersion objective. The 40x and the 100x objective are retractable (safeguarding against damages). With the 10x wide field eyepiece, you get a total magnification range 40x-1000x. The microscope has 360° rotatable monocular viewing and a coaxial coarse and fine focusing drive.

Included with your purchase:

DIN eyepieces (WF10x), DIN objectives (4x, 10x, 40x), mechanical desk with Nonius-scale, Abbe condenser, AC power adapter

Art.-Nr. · Art. No.:	51-02000
Modell · Model:	Erudit DLX
Typ · Type:	Durchlicht · Transmission
Vergrößerung · Magnification:	40x - 1000x
Okulare · Eyepieces:	DIN WF10x
Objektive · Objectives:	DIN 4x / 10x / 40x / 100x-Oil
Beleuchtung · Illumination:	LED Akku

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The Bresser TRINO RESEARCHER and BINO RESEARCHER are based on a solid metal body. The wide magnification range of 40x-1000x is reached with 4 objectives and one (changeable) pair of eyepieces. So, this instrument fulfils higher demands. The 40x- and the 100x-objective are retractable (safeguarding against damages). The 100x-objective is intended for oil immersion microscopy (immersion oil is included).

The halogen lamp is bright (20 W), and its intensity is variable with a dimmer.

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
Instead of a binocular head, the TRINO RESEARCHER has a trinocular head. With this equipment, the observer can make binocular observations and – at the same time! – take photographs or use a PC eyepiece.

Bino 57-22100

The adjustable Abbe condenser is equipped with a blue filter and an iris diaphragm. There is a focuser with coarse and fine drive for dead-sharp images, and a mechanical desk for exact object positioning. The microscope head is 360° rotatable. Individual adjustment of eye relief and diopter adjustment are also provided for.

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The new BIORIT ICD is the successful combination of good price, high quality and useful equipment. It is very suitable for viewing coins, stamps, minerals and many more. By the big working distance of 70mm and the wide focus range, also bigger objects can be viewed comfortable. An integrated illumination (12V / 10W) gives bright and even light onto the object.

All case parts are made of steel and aluminium; thus raises the Biorit ICD above the „toy“ class and give a solid base for the beginner and also for the experienced user.

Included with your purchase:

Stereo objective 2x, Pair of eyepieces WF10x, Spare bulb, Spare fuse, Dust cover

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The new ANALYTH ICD is a modern incident/transmitted light microscope with 360° rotatable microscope head and turnable zoom objective for different magnifications. The incident light source in the upper part of the microscope is for viewing opaque objects, the transmission light source is for translucent objects like preparation slides. It has two independently dimmable illuminations. By the big working distance of 70mm and the wide focus range, also bigger objects can be viewed comfortable.

All case parts are made of steel and aluminium; this gives a solid base for the for the experienced user.

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Included with your purchase:

Pair of eyepieces (WF10x), Stereo objectives (2x, 4x), Spare bulb, Spare fuse, Dust cover

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220 V / 240 V power plug at version with halogen illumination.

The Bresser RESEARCHER ICD is an all purpose instrument which is manufactured to the highest quality. It is a stereo microscope, but also has the option of being used as a biological microscope. This model has two sets of illuminators for each type of microscope style, which can be operated independently or together. The total intensity of the illuminators can be varied with a dimmer. The binocular head is 360° turnable and provides a comfortable view in all observing positions. The eyepieces can be set to the individual's eye relief and – with a diopter adjustment – to the user's eyesight. The microscope gives images, which are rich in details and in contrast.

The Halogen version (58-03000) offers powerful 10W illumination e.g. for cameras.

The LED version with integrated Accu is very versatile for outdoor use.

Included with your purchase:

360 degrees rotatable binocular head, diopter adjustment, 2 pairs of oculars, Illumination for transmitted and incident light, dimmer, glass- and black/white object table, 220 V / 240 V power plug.

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Eyepieces with ergonomic rubber eyecups

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The Bresser SCIENCE TRM-301 is a trinocular biological microscope of first-class quality: This means excellent imaging, good ergonomics, high stability and a modern attractive design. It is a laboratory microscope, which permits an intensive use, e.g. in medicine, biology, agriculture and industry. We can also recommend it to the enthusiast as his „microscope for a lifetime“.

A particularly important feature of the TRM-301 is the Köhler illumination. This optimal illumination – unbeatable bright, cool (going easy on the object), homogenous, reflex-free and rich in contrast – complements the superb quality of the microscope's optics. The Köhler illumination is standard for scientific research and microphotography. One gets it using the fixed field diaphragm and the height- and center adjustable Abbe condenser (numerical aperture 1.25; with iris diaphragm).

Included with you purchase:

lengthwise and cross moveable stage, ergonomic coaxial focusing knobs (concerning fine drive: 200 µm/full turn), continuous dimmer, plan-achromatic objectives (4x, 10x, 40x, 100x-oil retractable), 10x wide field eyepieces (mounted in 30° inclined eyepiece tubes), adjustment of eye relief, diopter adjustment and the phototube.

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
In the microscope's base, are a transmitted light halogen illuminator (20 W) as well as the power supply and continuous adjustable dimmers for reflected light and transmitted light. The trinocular head attachment works with both illumination modes. It also allows for photography and binocular visual observation. The inclination of the eyepiece tubes is 30°.

With the polarization head attachment, the instrument is suitable for the internal structure of materials and for condition monitoring of different mechanical apparatus in college and factory. The reflected light halogen illumination (20 W), which goes through the objective onto the object, can be modified for special purposes with a colour filter wheel and a polarizer/analyzer combination.

The Science ADL-601P represents excellent value for money.

Included with your purchase:

Wide field eyepieces (WF10x), planachromatic DIN objectives (4x, 10x, 40x retractable, 100x-oil retractable), Abbe condenser (numerical aperture 1,25; with iris diaphragm and auxiliary lens), immersion oil and a field diaphragm (in order to obtain Köhler illumination).

- 59-41981 DIN WF Okular 10x m. Fadenkr.
DIN WF eyepiece 10x w. crosshair
- 59-42100 Fotoadapter
Photo adapter

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The Science ADL-601F is the fluorescence version of the ADL-601P. In the microscope's base, there is a transmitted light halogen illuminator (20 W) as well as the power supply and continuous adjustable dimmers for reflected light and transmitted light. The trinocular head attachment functions with both illumination modes. It also allows for photography and binocular visual observation. The inclination of the eyepiece tubes is 30°.

The fluorescence attachment includes a high pressure mercury lamp (100 W) with external power supply, a filter wheel and a radiation shield. The mercury lamp produces the exciting radiation (namely ultraviolet radiation) in order to get fluorescence.

Included with your purchase:

Wide field eyepieces (WF10x), planachromatic DIN objectives (4x, 10x, 40x, 40x Glyc., 100x-oil retractable), Abbe condenser (numerical aperture 1,25; with iris diaphragm and auxiliary lens), immersion oil and a field diaphragm (in order to obtain Köhler illumination).

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
In contrast to normal light, polarized light waves vibrate in a single plane. Polarized light allows inspection of anisotropic objects, whose observable properties are dependant on the direction of the light waves, „in a new light“: crystals, animal hairs, feathers, muscles, nerve fibres and vegetable cell walls, all have a preferred direction of their biomolecules. These show structures when viewed using polarized light. Polarization microscopes are used particularly in mineralogy, geology, biology and medicine.


In order to change the object's direction in a controlled way, the Science MPO-401 has a round, center-adjustable and full rotatable object stage with degree scale. The illumination and observation unit encloses the halogen bulb lamp (20 W) with dimmer, the field diaphragm (used to achieve the Köhler illumination), the 360° rotatable polarizer (for the production of polarized light), centering implements for the objectives and the analyzer attachment (for the analysis of the polarized light). The analyzer attachment holds the 360° rotatable analyzer and optional available slips. The eyepieces are placed at the trinocular head with an inclination of 30°. They have eye relief adjustment as well as diopter adjustment.

Included with your purchase:

Wide field dividing eyepiece (WF10x), polarization DIN objectives (4x, 10x, 40x, 100x-Oil retractable), Abbe condenser with iris diaphragm and immersion oil.

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An the Science IVM 401, the usual direction of illumination and observation is inverted: One illuminates the sample from above and observes it from below. This type allows for the microscopy of sediments and transparent living organisms, which are found at the bottom of a liquid medium. For example, plankton or protozoa at the base of culture bottles, which should not turned over, can be investigated. So, this special microscope is suitable for public water departments, research institutes and universities, for medical treatment, animal husbandry, agriculture, limnology etc.

The Science IVM 401 has a large fixed object plate with moving stage system, a powerful halogen bulb illumination (30 W) with dimmer and field diaphragm, and a fourfold nosepiece. The trinocular head is turnable and has 30° inclined eyepiece tubes and a photo tube. There are eye relief adjustment and diopter adjustment, as well as coaxial coarse and fine focusing knobs.

Included with your purchase:

Wide field eyepieces (WF10x), planachromatic DIN objectives (10x, 25x, 40x), condenser (working distance 30 mm), filters (blue, green, ground glass) and specimen holders for rectangular or round culture utensils.

- ① Halogenbeleuchtung · *Halogen Illumination*
- ② Kondensator · *Condensor*
- ③ DIN-Objektive · *DIN Objectives*
- ④ DIN-Okulare · *DIN Eyepieces*
- ⑤ Fototubus · *Photo tube*
- ⑥ Dioptrieneinstellung · *Dioptre setting*
- ⑦ Objektstisch mit Objektführer
Object plate with stage moving system

- 59-41981 DIN WF Okular 10x m. Fadenkr.
DIN WF eyepiece 10x w. crosshair
- 59-42100 Fotoadapter
Photo adapter

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The Bresser SCIENCE MTL-201 is a special reflected light microscope, which allows the observation of non-transparent surfaces at higher magnifications (50x-800x). This solidly designed instrument is ideal for working in metallurgy (e. g. microscopy of polished samples of metals or alloys), mineralogy, precision engineering and electronics. It is suitable for scientific research, teaching demonstration and industrial production (e.g. quality controlling).

The MTL-201 is equipped with a large stage (250 mm x 230 mm), which can be moved coarsely by a quick adjustment lever and finely by a transverse and a longitudinal knob. The focuser has coarse and fine drive (coaxial; circa 70 µm/ full turn). The illumination is supplied by a halogen lamp (20 W) with a dimmer. The light comes through the objective from above onto the sample. The trinocular head has two 30° inclined ocular tubes and the photo tube. The photo tube allows for digital imaging, with an optionally available PC eyepiece. There is also a polarizer/analyzer combination and a field diaphragm, which can be centered; diopter and eye relief settings can be adjusted individually.


Included with your purchase:

Wide field eyepieces (WF10x), planachromatic DIN objectives (5x, 10x, 20x, 40x, 80x), filters (yellow, green, blue, ground glass).

- 59-41981 DIN WF Okular 10x m. Fadenkr.
DIN WF eyepiece 10x w. crosshair
- 59-42100 Fotoadapter
Photo adapter

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The Bresser SCIENCE ETD-101 is a high-quality stereo microscope with transmitted light function. It can be used for observation studies in medicine and health, in farming and forestry, in public security departments, schools and scientific research institutes as well as for inspection, assembling and repair of tiny parts in electronics and precision machine industries. This broad spectrum is due to the comprehensive possibilities of the ETD-101.

The broad furnished base with object stage and the extensive range of adjustable height allow the observation of larger objects also. The instrument has two illuminations (with 10 W halogen lamp each), which can be dimmed separately and can be used in three illumination modes: Transmitted light, reflected light and mixed light illumination. The incident angle of the reflected light illumination can be varied and fixed in a particularly user-friendly way. A pair of changeable wide field eyepieces (WF10x) and a stereo zoom objective (0,7x-4,5x) allow for the magnification range of 7x to 45x. The upright and true sided images are bright, true colored, rich in details and contrast. Naturally, the individual eye relief and the diopter adjustment can be set. The 360° turnable binocular head bears 45° inclined ocular tubes.

Included with your purchase:

Wide field eyepieces (WF10x), stereo zoom objective (0,7x - 4,5x)

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One of the most intelligent developments in the last few years is the MicroOcular. Connect your microscope with your PC and see what you would normally see trough your eyepiece on your monitor. The included software allows the user to easily record, save and play with his observations.*

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Included with your purchase (both models):

ArcSoft Photolmpression image editing software, adapter ring, USB connection cable

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
Note:

Bresser microscopes marked with this symbol can be used with a MicroOcular.

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Microscope ABC's

Aperture, numeric

Details of numerical aperture are engraved on the metal housing of every DIN lens (0.10, 0.25, 0.40 etc.).

These figures give the resolution capacity of the lens. Microscopic magnification should not exceed 1,000 times the aperture.

See: Resolution Capacity

Barlow Lens (BIOLUX STL)

The Barlow Lens, which is between the eyepiece and the monocular lens in the 51-15000, 51-16000 and 51-19000 models, enables you to adjust the magnification continuously. By pulling it out of the eyepiece tube, magnification is increased by more than half the value set.

Binoculars

A binocular sight enables you to view a slide with both eyes. This exploits the user's entire vision capacity.

Chemicals

Various chemicals and alcohol solutions are used in preparing specimens for lengthy storage. Xylol is used to brighten specimens.

Canada Balsam is a resin with the same optical refractive power as glass. After dewatering, staining and brightening a specimen the transparent slide is enclosed watertight in a drop of Canada balsam.

Slides thus treated are indefinitely preserved under normal circumstances. Distilled water serves to wash out the slides after staining or fixing.

See: Permanent Slide Cultures

Chemicals should be kept away from children.

Condenser

A condenser consists of one or more optical lenses. The lens bundles light to a point and thus ensures better illumination of slides.

Designs

There are two basic designs of microscope: direct and transmitted light. They account for by far the majority of

all microscopes used. The transparent slide is viewed using the transmitted light procedure. The maximum magnification here is limited to 3,000-4,000x. Non-transparent specimens (stones, coins) are viewed using the direct light procedure.

Direct-light microscopy

See: Designs and macro observation

Electron microscope

This type of microscope is mainly used in research and scientific laboratories.

In this type electron rays are bundled in a beam. The microscope image is made visible on a screen. Magnifications of up to 100,000x can be achieved with this type of microscope.

Eyepiece

The eyepiece is the lens system facing the eye in a microscope. Magnification of the lens used determines, together with that of the eyepiece, the overall magnification of the microscope.

The lifelike image generated by the lens is magnified by the eyepiece and viewed as a virtual image.

Focal length and focal point

The focal length is the gap between slide and eye. If it can be clearly seen it's called the focal point.

The normal space between the eyes when using the naked eye is 25 cm.

Magnifying glasses and microscopes have a shorter focal length (smaller gap) - the smaller it is, the greater the magnification.

Iris diaphragm

This is in the condenser and can be opened or closed with a lever. At higher magnifications the iris can be used to correct overexposure or too soft images.

Lens

The lens is the microscope lens system facing the slide. Lens and eyepiece magnification together determine microscope overall magnification. The lens as the „generative“ system renders a lifelike image. For users demanding higher standards and for photography our accessories range (pages 10 + 11) includes planachromatic and apochromatic lenses.

Microscope ABC's

Lens nosepiece

The rotatable nosepiece is used on microscopes with several lenses to switch between them and thus set varying magnifications.

Lighting

A light source is needed for both types. This source is either electrically (mains or battery power) or mirror generated.

Macro observation

Observing specimens such as coins, stamps or diamonds using the direct light procedure with a relatively low magnification (10x - 100x).

Magnification

Microscope magnification depends on that of the eyepiece and the lens.

Example: eyepiece 10x x lens 100x = magnification 1,000x

The useful magnification of a microscope lens is 1,000x that of the numerical aperture.

In economically priced student and hobby microscopes no data on lens aperture is given. The useful magnification in these types is approximately 300x400x - even if a final magnification of 900x or 1200x is feasible.

See: Aperture, Numerical

Mat and colour filters

These are inserted in the condenser filter holder and used to observe very bright or transparent slides.

Mechanical plate

This permits precise movement or fixing in place of a slide under the microscope.

Movement is by fine adjustment screws in two directions.

Metrology eyepiece

This is a removable DIN eyepiece with inbuilt 0-10 mm scale.

Micro eyepiece

Some models are supplied with these or can be fitted with them. This type of eyepiece allows you to connect microscope and PC and thus look at your slides on your monitor screen.

The micro eyepiece is inserted instead of the normal one and connected to the PC via USB. The software supplied allows you to easily make records of your observations and to process and store them.

These micro eyepieces fit in all 23 and 30 mm eyepiece holders.

Micro Projektor

Some models are supplied with a micro projector. The microscope image is projected on a screen or white wall and can then be viewed by several people at once.

Micrometer

The micrometer procedure in direct (OMO) and transmitted light (OMP) microscopy is used for precise evaluation and metrology down to the mm. In the process metrology is done by interaction between microscope eyepiece (adjustable) and microscope lens.

Microscope accessories

Microscope accessories include on the one hand those for the microscope itself - interchangeable eyepieces and lenses, condenser, camera adaptor, micro projector, etc.

On the other they include those for making slides - chemicals, staining media, slides themselves, glass plates, manual microtomes, microscopy aids, etc.

Microscope, optical system

A microscope has two lens systems - the eyepiece shows an enlarged image of the slide.

The lens shows a reversed image of it and this can be seen further enlarged through the eyepiece.

Monocular

A monocular insert lets you view the microscope image with just one eye.

Oil Immersion

In standard observations using higher magnification (90x/100x) it is usually impossible to improve image focus.

This lack is compensated for using oil immersion. A drop of immersion oil is added to the slide and the 90x or 100x lens then moved directly over it. The image is then seen through the oil and can easily be focussed.

Microscope ABC's

Permanent slide cultures

Permanent slide cultures are microscopically treated specimens glued to slides. The slide is sealed airtight and remains in unchanged condition for years.

Our professional permanent slide cultures are made in scientific laboratories and constantly monitored.

Photography (Biolam devices only)

To record successful microscopic observation same can be photographed using the focal or projection (with eyepiece) procedure. An SLR camera and suitable T2 ring plus a photo adaptor or accessory are needed for the purpose. The accessory or adaptor serves as connector between camera and microscope. See: T2 Ring

Plano concave mirror

A Plano-concave lens is a divergent or negative lens, thicker at the edge than in the centre and flat on one side.

It forms a virtual image and is frequently used to extend the effective focal length of a converging lens when they are used as a pair. It generates the light needed in transmitted light microscopy.

It uses an external light source (window, lamp) to guide a bundled light beam to the microscope. The light is bundled here by the parabolic mirror surface.

Resolution capacity

This is the ability to view two dots that are separate in the slide visually or view them separately at the photographic level. It is decisive to the question of what details can be seen at what magnification.

Stains

Microscope stains such as haematoxy-line, methylene blue, fuchsine, neutral red etc. serve to stain colourless slides. Staining makes certain details of a slide visible that would normally be colourless and hence invisible (e.g. starch granules).

Stereo Microscope

A stereomicroscope has two eyepieces making it possible to view the slide stereoscopically. Stereomicroscopes are mainly used for observing larger specimens not permeable to light (stones, leaves,

coins, etc.). The necessary space between lenses limits magnification to a maximum of approximately 100x.

T2-Ring

When photographing using a SLR camera the T2 ring connects the photo adaptor or accessory via an M42 thread with the camera.

Transmitted light microscopy

See: Designs

Trinocular

A trinocular accessory allows you to look into the microscope using both eyes. Photographs can be taken simultaneously.

Using the microscope

Two decisive errors are often made when using microscopes:

1. Too high a magnification is set.
2. The specimen is destroyed by incorrectly setting the lens.

Re 1.: For beginners a magnification between 50x and 300x suffices for simple transparent sectioned slides. Higher magnification is only needed when examining slides sliced with a microtome which are hence very thin. Very high magnifications are also used when examining blood specimens.

Re 2.: With high magnification, image clarity can only be adjusted shortly before the lens touches the slide. For correct adjustment the lens is therefore moved just above over the specimen. Image clarity is then carefully set whilst looking through the eyepiece.

Zoom eyepiece

The zoom eyepiece makes it possible to continuously adjust the focal length and hence magnification. Troublesome eyepiece change is unnecessary in microscopes with zoom lenses.

The focal length range with these eyepieces is 10x - 30x. The image quality of an eyepiece with fixed magnification is always better, though, than that of a zoom eyepiece. Zoom eyepieces are therefore usually used in student and hobby microscopes.

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Company Profile · Index

BRESSER was founded by Josef Bresser in 1957. The company began, by specializing in the import and distribution of binoculars.

In 1979, when Mr. Josef Bresser died, his son, Rolf Bresser, took over the business. He saw his main objective in the technical development and design of new and innovative optical products.

The product portfolio of BRESSER has been continuously expanded. In the last decades quite a number of patented product solutions have been established in all optical areas. This is why you can be sure that with the purchase of a BRESSER product you have purchased a state-of-the-art instrument at each price point.

Today BRESSER is one of the leading suppliers of optical instruments all over Europe and one of the largest in the world.

In order to remain a leader in the optical field, Rolf Bresser sold his business in 1999 to the American market leader MEADE INSTRUMENTS CORP. in Irvine, California. MEADE has established an excellent name in the telescope industry.

Founded in 1972 as a small mailorder company, MEADE – a publicly traded company on the NASDAQ stock exchange – is a world market leader in this market. Because of the advanced computerized MEADE Telescopes, new solutions were offered not only to amateur astronomers but also to novice astronomers. Technology such as this, had previously only been available to astronomy professionals.

In 2002, MEADE expanded by acquiring SIMMONS OUTDOOR CORP., which in-

cludes the Simmons and Weaver Brand names. This acquisition further established MEADE Instruments within the optical sports, outdoor and hunting businesses.

In the year 2003, the BRESSER product range was completely re-designed, expanded and re-engineered. Alongside the brands BRESSER and MEADE are new advanced programs from SIMMONS and WEAVER. The vast assortment allows anybody interested in optics to satisfy her or his personal requirement.

Also in 2003, the Josef-Bresser-Observatory was opened, equipped with the largest series-production telescope of the world, the MEADE LX200 16". This reference level observatory is not just open to MEADE customers and partners, but also to the public.

In 2005 MEADE expanded once again with the purchase of CORONADO, a world leader in ultra-narrowband filters used for sun observation.

Due to the continuous growth of MEADE and BRESSER the facilities at our location in Borken became too small. The European center moved to new modern premises to the neighbouring town of Rhede in 2006.

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
© 2008 Meade Instruments Europe GmbH & Co. KG.

No reproduction without permission.

All rights reserved. Subject to change. Errors excepted.

® The name „Bresser“ and the Bresser Logo are registered trademarks
of the Meade Instruments Europe GmbH & Co KG, Rhede, Germany.
The name „Meade“ and the Meade 4M Logo are registered trademarks
of the Meade Instruments Corporation, Irvine, USA.

KAT9936000DEGB0208BRESSER

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Warenzeichen der Meade Instruments Europe GmbH & Co KG,
sind eingetragene Warenzeichen der Meade Instruments Corporation,